



Product Information

Lecithin:cholesterol acyltransferase (LCAT), purified recombinant

Catalog No. RB-PRHL

Description	LCAT is an active protein present in normal human plasma involved in lipoprotein metabolism. This product is useful as a standardized source of LCAT for mass or activity assays.
Source	Recombinant protein from baby hamster kidney (BHK) cells
Quantity	~100 µg
LCAT Concentration	Lot specific
Storage and Handling	Store at -80° C. Avoid freeze/thaw cycles. Shelf life: 3 years, if stored properly. Wear appropriate gloves, protective clothing and eyewear. Follow safe laboratory practices.

Related Products

RB-LCAT	Roar LCAT Activity Assay Kit
RB-LIAG	LCAT Inhibitory Antiserum

For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

Roar Biomedical, Inc., Audubon Biomedical Center, 3960 Broadway, New York, NY 10032 USA
Tel: +1 (212) 280-2983 ▪ Fax: +1 (212) 280-2968 ▪ info@roarbiomedical.com ▪ www.roarbiomedical.com

©2004-2008 Roar Biomedical, Inc. All rights reserved. This information is subject to change without notice.